

LGX-50

LIQUID GAS EXTRACTOR


***The most effective way to automate
your liquid gas analysis***





ADVANTAGES

- No more manual prep
- Automatic internal standard addition – liquid or gas
- Sample is never exposed to the atmosphere
- System can also be used for VOC screening.

EST  **Rely On Us**
analytical 



Methane Curve



Gas Standard Overlay

APPLICATION RESULTS

The saturated gas curves were linear with results of 0.999 or better linear regression. Carryover after the high standard was less than 0.10% and the precision of the IS and the gases tested was better than 10% RSD.

OVERALL VALUE

- Best Value on the Market
- Proudly made in the USA

BUY WITH CONFIDENCE

- EST analytical is a worldwide leader in GC sample introduction

EST  **Rely On Us**
 analytical 